


Tristel Protect™

Transportation bags for instruments
with integrated traceability



Transportation bags for instruments with integrated traceability

Tristel Protect bags are designed specifically for the transportation of instruments, with traceability in mind. Each bag is identified with a unique and patented label, incorporating a resealable pocket. This pocket enables you to store traceability data, without compromising the integrity of the instrument inside the bag.



Why choose Tristel Protect?

A tamper-proof system

The Tristel Protect pocket label is patented and comprises a special seal material. When the instrument inside the bag is clean, the seal is intact (solid blue). Upon opening the bag, the seal transfers onto the bag to show an "OPEN VOID" marking. This indicates that the bag can only be used to transport a dirty instrument and minimises the risk of mistakes in handling instruments.

No capital investment

Tristel Protect requires minimal capital investment when compared to a traditional tray transportation system.

Two sizes to accommodate a wide range of instruments

Large Tristel Protect bags each measure 30 x 70 cm. They can be used for the transportation of:

- Nasendoscopes
- Intubation endoscopes
- Transvaginal ultrasound probes
- Transrectal ultrasound probes

Small Tristel Protect bags each measure 23 x 33 cm and can be used for the transportation of:

- Laryngoscope blades
- Manometry catheters
- Pachymeters
- Tonometer prisms

Small footprint

Tristel Protect requires minimal space. Each box only measures 29 x 29 x 12cm.

A strong product

Tristel Protect bags are made of 100% Low Density Poly Ethylene to prevent rips and tears.

How to use Tristel Protect



1 Take a Tristel Protect bag from the box.



2 Place the clean instrument inside the bag and seal it. The seal is intact (solid blue) to confirm that the instrument inside the bag is clean.



3 Press and peel the corner of the label to open the pocket.



4 Secure traceability data within the marked area and reseal the pocket. The instrument is now ready to be transported to the patient.


Upon arrival to the patient, remove traceability data from the pocket and attach it to patient records.



5 Open the bag as shown. The seal transfers onto the bag to show an "OPEN VOID" marking. Remove the instrument from the bag and use it on the patient. Do not throw the bag away.

Re-use the Protect bag to transport the dirty instrument back to the decontamination room. The pocket can be used again to complete the traceability between device, disinfection procedure and patient. Upon arrival to the decontamination room, the used Protect bag is discarded to clinical waste.

The next logical step in the decontamination cycle



Tristel Protect represents the logical next step following decontamination with the Tristel Trio Wipes System. A system comprising cleaning, high-level disinfection, rinsing and transport can be established. At the heart of this system is traceability.



The Protect carton is designed in such a way, that the Trio carton can be placed on top of it to create a complete system.

For more information about the Tristel Trio Wipes System please visit www.tristel.com/tristel-products/tristel-trio-wipes-system.

Product options

TSLo71001 Small Tristel Protect Bags

1 box of 50 small bags



TSLo70901 Large Tristel Protect Bags

2 boxes of 50 large bags



Tristel Protect patent pending:
United Kingdom: GB 2526313 & International: WO 2015/177518

Tristel

Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

Australia: Tristel Pty Ltd, 144 Church Street, Brighton, Victoria, 3186
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

New Zealand: Tristel New Zealand Limited, 23 Birch Avenue, Judea, Tauranga P.O. Box 3110
T +64 (0)7 5771560 - F +64 (0)7 5771567 - E info@tristel.co.nz

Hong Kong & Taiwan: Tristel Asia Limited, 2404 Worldwide House, 19 Des Voeux Road Central, Hong Kong
T +85 (0)23667 9253 - F +85 (0)23667 9002 - E customerservice@tristel.com.hk

Copyright © Tristel Solutions
Mkt-Bro-098-2
May 2017

