

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model
Alpinion Medical	Skin Surface Transducers	No Series Specified	MP1-5X, SP1-5, SP1-5i, SP1-5T, SP1-5X, SP3-8, SP3-8T, SC1-4H, SC1-4HS, SC1-6, SC1-6H, C1-6, C1-6T, C1-6i, C1-6i, C5-8, C5-8N, C5-8NT, L3-12, L3-12H, L3-12HWD, L3-12i, L3-12T, L3-12X, L3-8, L3-8i, L8-17, L8-17i, L8-17X, CW2.0, CW5.0
	Endocavity Transducers		E3-10, E3-10H, EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X
	Intraoperative Transducers		IO3-12, IO8-17, IO8-17T
BK Medical / BK Ultrasound	Skin Surface Transducers	bk3000/bk5000	5P1, 6C2, 6C2s, 9C2, N13C5
	Endocavity Transducers		8L2, 18L5, 14L3, 13L4W, 10L2W
	Intraoperative Transducers		E10C4
	Intraoperative Transducers		E14C4t, E14C4, E14CL4b, 20R3
	Endocavity Transducers	Flex Focus & Pro Focus	l14C5T, N11C5s, l12C4f, X18L5s, l14C5l, l12C5b, X12C4
	Endocavity Transducers		8862, 8863, 8836, 8814, 8815, 8816, 8824
	Ultrasound Unit	No Models Specified	2050, 2052, 8808, 8808e, 8818, 8838, 8848
Carestream	Skin Surface Transducers	No Series Specified	bk3000
	Endocavity Transducers		bk3500
Esaote	Endocavity Transducers	No Series Specified	bk5000
	Intraoperative Transducers		Flex Focus
	Veterinary Endocavity Transducer		No Models Specified
GE Healthcare	Skin Surface Transducers	Invenia ABUS	6C2, 6C2s, 9C2, 14L3, 13L4W, 10L2W, 18L5, N13C5, 5P1
	Endocavity Transducers	Logiq Systems - E9	E14C4t, E14C4, E14CL4b, E10C4, 20R3
	Intraoperative Transducers		EC123, EC1123
	Skin Surface Transducers	Logiq Systems - S7	IOE323
	Endocavity Transducers		LA vet probe
	Endocavity Transducers	No Models Specified	C1-5-D, C1-6-D, C2-9-D, M6C-D, 9L-D, 11L-D, ML6-15-D, 6S-D, RAB2-5-D, RAB4-8-D
Endocavity Transducers	No Models Specified	IC5-9-D, RIC5-9-D	
Endocavity Transducers	No Models Specified	L8-18i-D	
Endocavity Transducers	No Models Specified	C1-5-D, 8C, 9L-D, 11L-D, ML6-15, L8-18i-D, RAB4-8-D	
Endocavity Transducers	No Models Specified	IC5-9-D	

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model	
GE Healthcare	Skin Surface Transducers	Logiq P Series	4C, 5CS, 8C, 8L, 9L, 11L, ML6-15, 5S, 6S, 7S, 4D3C-L, 4D5C-L	
	Endocavity Transducers		BE9CS, E8C, E8CS, 4DE7C	
	Intraoperative Transducers		I739, T739	
	Skin Surface Transducers	Logiq C Series	C2-5-RC, L6-12-RC, 3S-RC	
	Endocavity Transducers		E8CS-RC	
	Skin Surface Transducers	Logiq A Series	3.5CS, 5CS, 8C, 8L, 9L, 10L, 11L, 12L, 3S, 5S, 7S	
	Endocavity Transducers		BE9CS, E8CS	
	Intraoperative Transducers		I739, T739	
	Skin Surface Transducers	Logiq E Series	4C-RS, 8L-RS, 9L-RS, 12L-RS, 3S-RS, 6S-RS	
	Endocavity Transducers		E8C-RS	
	Intraoperative Transducers		8C-RS, I12L-RS, I739-RS, T739-RS	
	Skin Surface Transducers	Vivid E9	4C-D, 9L-D, 11L-D, ML6-15-D, 6S-D, 12S-D, 4VD	
	Skin Surface Transducers	Vivid P Series	4C, 5CS, 8C, E8CS, 8L, 11, T739, 3S, 6S	
	Endocavity Transducers		E8CS	
	Intraoperative Transducers		T739	
	Skin Surface Transducers	Vivid S/e/i/q Series	4C-RS, 8L-RS, 9L-RS, 12L-RS, 3S-RS, 3Sc-RS, 5S-RS, 6S-RS, 7S-RS, 12S-RS, M4S-RS	
	Endocavity Transducers		E8C-RS	
	Intraoperative Transducers		8C-RS	
	Skin Surface Transducers	Voluson Expert Series	4C-RS, 8L-RS, 9L-RS, 12L-RS, 3S-RS, 3Sc-RS, 5S-RS, 6S-RS, 7S-RS, 12S-RS, M4S-RS	
	Endocavity Transducers		E8C-RS	
	Intraoperative Transducers		8C-RS	
	Skin Surface Transducers	Voluson Performance Series	4C-RS, 12L-RS, 3Sc-RS, RAB2-5-RS	
	Endocavity Transducers		E8C-RS, RIC5-9W-RS	
	Skin Surface Transducers	Voluson Signature Series	4C-RS, AB2-7-RS, C1-5-RS, 9L-RS, 12L-RS, ML6-15-RS, 3Sc-RS, RAB2-5-RS, RAB4-8-RS,	
	Endocavity Transducers		E8C-RS, RIC5-9W-RS	
	Intraoperative Transducers		8C-RS	
	Skin Surface Transducers	Venue 40 Series	4C-SC, 12L-SC, L8-18i-SC, 3S-SC	
	Endocavity Transducers		E8CS-SC	
	Healcerion	Mobile Ultrasound Device	No Series Specified	SONON 300C

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model	
Hitachi	Ultrasound Unit	Arietta 60	No Models Specified	
		Arietta 70		
		F31		
		F37		
		Hi-Vision Ascendus		
		Hi-Vision Avius		
		Hi-Vision Preirus		
		Noblus		
		ProSound Alpha 6		
		ProSound Alpha 7		
	ProSound F75			
	Ultrasound Transducers	No Series Specified	C22I, C22K, C22P, C22T, C25, C251, C252, C25P, C35, C41, C41B, C41L47RP, C41RP, C41V, C41V1, C42, C42K, C42T, CC41R, CC41R1, CL4416R, L34, L43K, L44, L441, L44K, L44LA1, L46K, L46K1, L51K, L53K, L53K1, L55, L64, S12, S121, S21, S211, S22, S31, S31KP, S42, SML44, VC41V,	
			ASU-1000B-3.5, ASU-1000C-3.5, ASU-1001, ASU-1002, ASU-1003, ASU-1005, ASU-1009, ASU-1010, ASU-1011, ASU-1012, ASU-1013, ASU-35-3, ASU-35CWD-2, ASU-35WL-10, ASU-35WSJ-7.5, ASU-36WL-10, ASU-64B, ASU-65B, ASU-67	
EUP-B314, EUP-B512, EUP-B514 EUP-B712, EUP-B715, EUP-C321T, EUP-C511, EUP-C514, EUP-C516, EUP-C524, EUP-C532, EUP-C715, EUP-CV714, EUP-F334, EUP-F531, EUP-L34T, EUP-L52, EUP-L53, EUP-L53L, EUP-L53S, EUP-L54MA, EUP-L73S, EUP-L74M, EUP-L75, EUP-LS85, EUP-O33JA, EUP-O53T, EUP-O54I, EUP-O732T, EUP-S50A, EUP-S52, EUP-S63, EUP-S70, EUP-S72, EUP-S80				

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Ultrasound Transducers	No Series Specified	UST-2265-2, UST-2266-5, UST-5024N-3.5, UST-5037P-3.5, UST-5039P-3.5, UST-5043P-3.5, UST-5044-3.5, UST-5045P-3.5, UST-5046-3.5, UST-5047-5, UST-52103, UST-52105, UST-52108, UST-52109, UST-52114P, UST-52124, UST-52127, UST-52129, UST-5226D-3.5/2, UST-5236D-2.5/2, UST-5242D-2.5/2, UST-5244-3.5, UST-5255-2.5, UST-5256-3.5, UST-5257-5, UST-5263D-2.5/2, UST-5265-2.5, UST-5266-3.5, UST-5267-5, UST-5268P-5, UST-5274-2.5, UST-5275-3.5, UST-5276-5, UST-5279-3.5, UST-5281-5, UST-5283-2.5, UST-5284-2.5, UST-5285-3.5, UST-5286-2.5, UST-5287-3.5, UST-5292-5, UST-5294-5, UST-5295, UST-5296, UST-5297, UST-5298, UST-5299, UST-533, UST-534, UST-535, UST-536, UST-5410, UST-5410H, UST-5411, UST-5412, UST-5413, UST-5415, UST-5417, UST-5419, UST-547, UST-5512U-7.5, UST-5514DTU-7.5, UST-5515-7.5, UST-5521L-7.5, UST-5523L-7.5, UST-5524-5, UST-5524-7.5, UST-5526L-7.5, UST-5531, UST-5534T-7.5, UST-5539-7.5, UST-5540P-7.5, UST-5542, UST-5543, UST-5545, UST-5546, UST-5548, UST-5551, UST-5561-7.5, UST-556T-7.5, UST-556TU-7.5, UST-567, UST-568, UST-5710-7.5, UST-5711-7.5, UST-5712, UST-5713T, UST-574T-7.5, UST-575-7.5, UST-576T-7.5, UST-577-7.5, UST-579T-7.5, UST-5818-5, UST-5819T-5, UST-657-5, UST-660-7.5, UST-665P-5, UST-667P-5, UST-670P-5, UST-675P, UST-676P, UST-677P, UST-9101-7.5, UST-9102-3.5, UST-9102U-3.5, UST-9103-5, UST-9104-5, UST-9111-5, UST-9112-5, UST-9113P-3.5, UST-9114-3.5, UST-9115-5, UST-9116P-5, UST-9118, UST-9119, UST-9120, UST-9121, UST-9123, UST-9124, UST-9125-7.5, UST-9126, UST-9127, UST-9128, UST-9130, UST-9132I, UST-9132T, UST-9133, UST-9135P, UST-9136, UST-9136U, UST-9137, UST-9145, UST-9146I, UST-9146T, UST-9147, UST-9148, UST-9149, UST-9151, UST-932-5, UST-932-7.5, UST-935N-5, UST-939-3.5m, UST-939D-3.5, UST-943-3.5, UST-944B-3.5, UST-945B-5, UST-945BP-5, UST-952-5, UST-953P-5, UST-954-5, UST-957-7.5, UST-961-5, UST-961-7.5, UST-964P-5, UST-968P-5, UST-969-5, UST-970-7.5, UST-972- 3.5, UST-974-5, UST-975-7.5, UST-977-5, UST- 978-3.5, UST-979-3.5, UST-980-5, UST-980P-5, UST-981-5, UST-981P-5, UST-984-5, UST-984P-5, UST-985-7.5, UST-987-7.5, UST-988-2.5, UST- 989-3.5, UST-989D-3.5, UST-990-5, UST-992-5, UST-994P-5, UST-995-7.5, UST-996P-5, UST-MC11-8731
Philips	Endocavity	Affiniti	3D9-3v
	Skin surface		C5-1
			C6-2 Compact
			C8-5 Compact
Intraoperative	L12-3		
L12-4 Compact			
L12-5 Compact			
L15-7io			

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model
Philips	Skin surface	Affiniti	S12-4 Compact
			S4-2 Compact
			S5-1
			S8-3 Compact
			V6-2
	Endocavity	ClearVue	3D9-3V
	Skin surface		C-2
			L12-4
			S4-1
			V6-2
	Endocavity	EPIQ	3D9-3V
	Skin surface		C5-1
			C8-5 Compact
			L12-3
			L12-5 Compact
	Intraoperative		L15-7io
	Skin surface		S12-4 Compact
			S5-1
			S8-3 Compact
			V6-2
			X5-1
	Endocavity	HD11	3D9-3V
			C5-2
			L12-3
	Intraoperative	L15-7io	
	Skin surface	HD11 XE	V6-2
	Endocavity	HD15	3D9-3V
	Skin surface		C5-1
			C5-2
			L12-3
	Intraoperative		L15-7io
	Skin surface	V6-2	
Endocavity	iU22	3D9-3V	
Skin surface		C5-1	
		C5-2	

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model
Philips	Intraoperative	iU22	L15-7io
	Skin surface		S4-1
			S5-1
			V6-2
			X7-2
		CX50	C5-1
	C8-5 Compact		
	L12-3		
	L12-5 Compact		
	L15-7io		
	S12-4 Compact		
	S5-1		
	S8-3 Compact		
	iE33	C5-1	
		C5-2	
		L15-7io	
		S5-1	
	Skin surface	X5-1	
		X7-2	
		Sparq	C5-1
			C6-2 Compact
	L12-4 Compact		
	L15-7io		
	Skin surface	HD5	S4-2 Compact
			C5-2
			L12-3
			S8-3 Compact
HD6		V6-2	
		C5-2	
		L12-3	
		L12-3	
Lumify	HD7	C5-2	
		L12-4	
		S4-1	
		Xperius	Intraoperative

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model	
Philips	Skin surface	Xperius	C5-2	
		CX30	C8-5 Compact	
	L12-4 Compact			
	L15-7io			
	S4-2 Compact			
	S8-3 Compact			
	Intraoperative	No series specified	3D6-2	
			3D8-4	
	7,5MHZ Endo			
	BP10-5ec			
	C10-3v			
	C10-4ec			
	C6-3			
	C8-4v			
	C8-5			
	C9-2			
	Skin surface	No series specified	C9-3io	
			C9-3v	
	C9-4			
	C9-4ec			
	C9-4v			
	C9-4v Compact			
	C9-5ec			
	Endocavity		No series specified	L11-3
				L12-5
				L17-5
		L18-5		
		L8-4		
		L9-3		
		S12		
S12-4				
S3-1				
S4-2				
S5-2				
S8				
S8-3				
V8-4				
VL13-5				
Skin surface	No series specified	X3-1		
		X6-1		

TRISTEL DUO ULT MATERIAL COMPATIBILITY WITH ULTRASOUND MEDICAL DEVICES

Manufacturer	Instrument Type	Instrument Series	Model
Samsung	Endocavity Transducers	No Series Specified	ER4-9/10ED
			V5-9
			3D4-9(ES)
			EVN4-9/10ED
	Skin Surface Transducers	No Series Specified	CA1-7A
			C2-6(IC)
			CA2-8A
			CF2-8
			CA2-8AD
			L5-13(IS)
			L7-16(IS)
			LA3-16AI
			L3-12A
			LA3-14AD
L5-13/50			
Samsung	Skin Surface Transducers	No Series Specified	CV1-8A, CV1-8AD, C2-8(ED), CA2-9A, CA2-9AD, CA1-7AD, CA2-6BM, L3-8(ED), LA3-14AD, LA3-14A, PA1-5A, PA1-5AD, SP3-8
	Endocavity Transducers		3DC-6(EN), 3D4-9(ES), EV3-10B, EVN4-9, EC4-9(CD), VR5-9, E3-12A, EA2-11B
	Intraoperative Transducers		LS6-15
Supersonic Imagine	Endocavity Transducers	Aixplorer	SE12-3
	Skin Surface Transducers	Aixplorer	SL15-4
		Aixplorer	SC6-1
Verathon	Skin Surface Transducers	BladderScan	BladderScan Prime Plus

Disclaimer: Material compatibility testing performed via different methods including wiping and soaking. Tested materials subjected to varying number of cycles.

For more information contact Tristel team on cooperation@tristel.com



Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Australian sponsor: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

New Zealand: Tristel New Zealand Limited, P.O Box 15451, Tauranga, New Zealand
T +64 (0)7 5771560 - F +64 (0)7 5771567 - E info@tristel.co.nz

Hong Kong & Taiwan: Tristel Asia Limited, 21st Floor, 168 Electric Road, North Point, Hong Kong
T +852 2895 6968 - F +852 2869 4388

Copyright © Tristel Solutions
Mkt-Val-457-1
March 2018