

Tristel Fuse™

for Stella

COMPATIBILITY 

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY			
Manufacturer	Instrument Type	Instrument Series	Model
BK Ultrasound	Intraoperative transducer	Flex Focus & Pro Focus	8862
			8809
			8863
			8836
			8814
			8815
			8816
			8824
	Endocavity transducer		8826
			2050
			2052
			8808
			8808e
			8818
			8838
			8848
Carestream	Skin surface transducer	Touch Prime	6C2
			6C2s
			9C2
			N13C5
			5P1
	Intraoperative transducer		14L3
			13L4w
			10L2w
			18L5
			18L5s
	Endocavity transducer		E14C4t
			E14C4
			E14CL4b
			E10C4
	Ultrasound transducer		20R3

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Esaote	Intraoperative transducer	Not specified	IOE323
	TOE probe		TEE022
			TEE122
			ST2612
	ST2613		
Endocavity transducer	LA vet probe		
Fujifilm Fujinon	Endocavity transducer	Series 100	All models
		Series 200	
		Series 250	
		Series 310	
		Series 410	
		Series 450	
		Series 485	
GE Healthcare	TOE transducer	Not specified	TEE 6T
			TEE 6T-RS
			TEE 6Tc
			TEE 6Tc-RS
			TEE 6Tv
			TEE 9T
			TEE 9T-RS
Gyrus / ACMI / Olympus	Flexible cystoscope	ICN endoscope	Not specified
Hitachi	Intraoperative transducer	ARIETTA	C22I
			C22K
			C22P
			C22T
			C25P
			C42K
			C42T
			L43K
			L44K
			L44LA1
			L46K
			L46K1
			L51K
			L53K
			S31KP
UST-533			

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Skin surface transducer	ARIETTA	C251
			C252
			C35
			C41
			C42
			L34
			L44
			L441
			L55
			L64
			S12
			S21
			S211
			S22
			S31
	S42		
	Endocavity transducer	C41B	
		C41L47RP	
		C41RP	
		C41V1	
		CC41R	
		CC41R1	
		CL4416R	
		R41R	
		R41RL	
		VC41V	
	Intraoperative transducer	Noblus	EUP-054J
Skin surface transducer	S21		
Endocavity transducer	C25		
Skin surface transducer	LISENDO 88o	C41V	
Skin surface transducer		S121	

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Skin surface transducer	ProSound	ASU-1010
			ASU-1012
			ASU-1013
			UST-52105
			UST-52124
			UST-52127
			UST-52129
			UST-5296
			UST-5298
			UST-5299
			UST-5411
			UST-5412
			UST-5413
			UST-5415
			UST-5417
			UST-5419
			UST-5542
			UST-5543
			UST-5548
			UST-568
			UST-5712
			UST-9102U-3.5
			UST-9115-5
			UST-9123
			UST-9127
			UST-9128
			UST-9130
			UST-9133
			UST-9136U
			UST-9147
UST-987-7.5			
UST-996P-5			
CW Doppler	UST-2265-2		
UST-2266-5			

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model	
Hitachi	Intraoperative transducer	ProSound	UST-5045P-3.5	
			UST-52109	
			UST-52114P	
			UST-533	
			UST-534	
			UST-536	
			UST-547	
			UST-5534T-7.5	
			UST-5710-7.5	
			UST-5713T	
			UST-579T-7.5	
			UST-9104-5	
			UST-9120	
			UST-9132I	
			UST-9132T	
			UST-9135P	
			UST-9146I	
			UST-995-7.5	
			UST-MC11-8731	
			Endocavity transducer	UST-670P-5
	UST-672-5/7.5			
	UST-675P			
	UST-676P			
	Endocavity transducer		UST-677P	
			UST-678	
			UST-9118	
			UST-9124	
	Intraoperative transducer		HI VISION	UST-984-5
				EUP-B512
				EUP-B514
EUP-B712				
EUP-B715				
EUP-F334				
EUP-F531				
EUP-O53T				
EUP-O54I				
EUP-O732T				

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model		
Hitachi	Skin surface transducer	HI VISION	EUP-C511		
			EUP-C514		
			EUP-C516		
			EUP-C524		
			EUP-C532		
			EUP-C715		
			EUP-CV714		
			EUP-L34T		
			EUP-L52		
			EUP-L53		
			EUP-L53L		
			EUP-L53S		
			EUP-L54MA		
			EUP-L73S		
			EUP-L74M		
			EUP-L75		
			EUP-S70		
			EUP-S72		
			EUP-S80		
			Endocavity transducer	EUP-CC531	
	EUP-CC531S				
	EUP-R54AW-19				
	Endocavity transducer		EUP-R54AW-33		
			EUP-U533		
			EUP-V53W		
			EUP-V73W		
					EUP-VV731

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Ultrasound transducer	Not specified	L53K1
			SML44
			ASU-1000B-3.5
			ASU-1000C-3.5
			ASU-1001
			ASU-1002
			ASU-1003
			ASU-1005
			ASU-1009
			ASU-1011
			ASU-35-3
			ASU-35CWD-2
			ASU-35WL-10
			ASU-35WSJ-7.5
			ASU-36WL-10
			ASU-64B
			ASU-65B
			ASU-67
			UST-5024N-3.5
			UST-5037P-3.5
			UST-5039P-3.5
			UST-5043P-3.5
			UST-5044-3.5
			UST-5046-3.5
			UST-5047-5
			UST-52103
			UST-52108
			UST-5226D-3.5/2
			UST-5236D-2.5/2
			UST-5242D-2.5/2
			UST-5244-3.5
			UST-5255-2.5
			UST-5256-3.5
UST-5257-5			
UST-5263D-2.5/2			
UST-5265-2.5			
UST-5266-3.5			

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Ultrasound transducer	Not specified	UST-5267-5
			UST-5268P-5
			UST-5274-2.5
			UST-5275-3.5
			UST-5276-5
			UST-5279-3.5
			UST-5281-5
			UST-5283-2.5
			UST-5284-2.5
			UST-5285-3.5
			UST-5286-2.5
			UST-5287-3.5
			UST-5292-5
			UST-5294-5
			UST-5295
			UST-5297
			UST-535
			UST-5410
			UST-5410H
			UST-5512U-7.5
			UST-5514DTU-7.5
			UST-5515-7.5
			UST-5521L-7.5
			UST-5523L-7.5
			UST-5524-5
			UST-5524-7.5
			UST-5526L-7.5
			UST-5531
			UST-5539-7.5
			UST-5540P-7.5
			UST-5545
UST-5546			
UST-5551			
UST-5561-7.5			
UST-556T-7.5			
UST-556TU-7.5			
UST-567			

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Ultrasound transducer	Not specified	UST-5711-7.5
			UST-574T-7.5
			UST-575-7.5
			UST-576T-7.5
			UST-577-7.5
			UST-5818-5
			UST-5819T-5
			UST-657-5
			UST-660-7.5
			UST-665P-5
			UST-667P-5
			UST-9101-7.5
			UST-9102-3.5
			UST-9103-5
			UST-9111-5
			UST-9112-5
			UST-9113P-3.5
			UST-9114-3.5
			UST-9116P-5
			UST-9119
			UST-9121
			UST-9125-7.5
			UST-9126
			UST-9136
			UST-9137
			UST-9145
			UST-9146T
			UST-9148
			UST-9149
			UST-9151
			UST-932-5
			UST-932-7.5
			UST-935N-5
UST-939-3.5			
UST-939D-3.5			
UST-943-3.5			
UST-944B-3.5			

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Ultrasound transducer	Not specified	UST-945B-5
			UST-945BP-5
			UST-952-5
			UST-953P-5
			UST-954-5
			UST-957-7.5
			UST-961-5
			UST-961-7.5
			UST-964P-5
			UST-968P-5
			UST-969-5
			UST-970-7.5
			UST-972-3.5
			UST-974-5
			UST-975-7.5
			UST-977-5
			UST-978-3.5
			UST-979-3.5
			UST-980-5
			UST-980P-5
			UST-981-5
			UST-981P-5
			UST-984P-5
			UST-985-7.5
			UST-988-2.5
			UST-989-3.5
			UST-989D-3.5
			UST-990-5
			UST-992-5
			UST-994P-5
Intraoperative transducer		EUP-LS85	
		EUP-B314	
		EUP-O33JA	

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Hitachi	Skin surface transducer	Not specified	EUP-C321T
			EUP-S50A
			EUP-S52
			EUP-S63
	Endocavity transducer		EUP-U531
			UST-664-5/7.5
			UST-666-5/7.5
			UST-669-5/7.5
Karl Storz	Hysteroscope	Hysteroscope	KS00337
Philips	Endocavity transducer	iU22	C8-4v (p/n 4535611 7292x where x is 3 or greater and p/n 4535612 8750x where x is 2 or greater)
		HD15	
		HD11 XE	
		HD9	
		HD7 XE	
		HD7	
		HD6	
		HD5	
		CX50	
		Sparq	
		CX30	
		HD15	BP10-5ec
		HD11 XE	
		Affinity 50	
		Affinity 70	
		EPIQ7	C9-4v
		EPIQ5	
		Affinity 50	
		Affinity 70	
		CX50	
		EPIQ7	
		EPIQ5	
		iU22	
		Affinity 70	
		CX50	

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model	
Philips	TOE probe	iE33	S8-3t	
		EPIQ7/EPIQ5/Affinity70/ Affinity50/iE33/HD11 XE/HD11	S7-3t	
		iE33	T6207	
		iE33	D2cwc	
		EPIQ7		
		EPIQ5		
		iU22		
		CX50		
		HD11		
		HD9		
		HD15		
		TOE probe		HD11 XE
			X8-2t	
	X8-2ti			
	EPIQ7		D2tcd	
	EPIQ5			
	iE33			
	iU22			
	HD11 XE			
	EPIQ7		S8-3	
	EPIQ5			
	Affinity50			
	Affinity70			
	CX30			
	CX50			
	iE33			
	HD11			
HD11 XE				
HD15				

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Philips	TOE probe	iE33	S7-2t
		iU22	
		HD15	
		HD11 XE	
		HD11	
		HD9	
		HD7 XE	
		HD7	
		HD6	
		HD5	
	TOE probe	iE33	S7-2 Omni (21778A)
		iU22	
		HD15	
		HD11 XE	
		HD11	
		HD9	
		HD7 XE	
		HD7	
		HD6	
		HD5	
	TOE probe	HD5	S7-2 Omni (21778A)
		iE33	T6H/21378A OMNI III
		HD11 XE	
		HD7 XE	
		EPIQ7	X7-2t
		EPIQ5	
		Affinity70	
		iE33	
		CX50	
		Sparq	
Skin surface transducer	EPIQ7	S5-1	
	EPIQ5		
	Affinity70		
	CX50		
	iU22		
	iE33		

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Philips	Skin surface transducer	EPIQ7	S12-4
		EPIQ5	
		Affinity70	
		Affinity50	
		iE33	
		HD15	
		HD11	
		CX50	
		HD11 XE	
		iE33/iU22	
		EPIQ7/iE33	X5-1
		EPIQ7/iE33/iU22	X7-2
		Ultrasound transducer	Not specified
	Skin surface transducer	EPIQ7	C5-1
		QPIQ5	
		Affinity 70	
		iU22	
		iE33	
		HD15	
		iE33	C5-2
		iU22	
	Skin surface transducer	HD15	C5-2
		HD11 XE	
		HD11	
		HD9	
		HD7 XE	
		HD7	
HD6			
HD5			
Endocavity transducer	iE33	C9-4v	
	EPIQ7		
	EPIQ5		
	Affinity50		
	Affinity70		
	CX50		

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Philips	Skin surface transducer	EPIQ7	C8-5
		EPIQ5	
		Affinity 70	
		Affinity 50	
		iE33	
		iU22	
		HD7 XE	
		HD11 XE	
		HD9	
		HD11	
		HD15	
		CX30	
		CX50	
		Intraoperative transducer	
	iU22		
	HD11 XE		
	iE33		L9-3
	iU22		
	HD15		
	HD9		
	HD11 XE		
	iE33		L11-3
	EPIQ7		L15-710
	Intraoperative transducer	EPIQ5	L15-710
		Affinity 70	
		Affinity 50	
		iE33	
		iU22	
		HD11	
		HD15	
CX30			
CX50			

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model	
Toshiba	Accessories	Xario	PVU-781VT	
			PVU-674BT	
			PLU-1005BT	
			PLU-1204BT	
	Endocavity transducer	Nemio & Famio	PVG-366M	
			PVG-381M	
			PVG-681S	
			PVG-600S	
			PVG-601V	
			PVG-630V	
			PVM-651VT	
			PVM-656AT	
			PVM-740RT	
			PVL-715R	
			PVL-715RT	
			Aplio & Viamo	PVT-661VT
				PVT-674BT
				PVT-781VT
		PVT-781VTE		
		Xario	PVL-715RST	
			PVT-770RT	
			PVL-715RS	
	Intraoperative transducer		Aplio & Viamo	PVT-745BTF
				PVT-745BTH
				PLT-308P
				PLT-705BTF
		PLT-705BTH		
		PVT-745BTV		
	Intraoperative transducer	Aplio & Viamo	PVT-350BTP	
			PVT375ST	
PVT-375MV				
PVT-382MV				
PVT-675MV				
PVT-675MVL				
PVT-681MV				
PVT-681MVL				

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model		
Toshiba	Intraoperative transducer	Nemio & Famio	PVF-738F		
			PVF-738H		
			PVF-745V		
			PLF-308P		
			PLG-308P		
			PLM-1202S		
			PLM-1204AT		
		Xario	PVU-745BTV		
			PVU-745BTF		
			PVU-745BTH		
			PVU-770ST		
			Skin surface transducer	Xario	PVU-681MVL
			Pencil transducer	Aplio & Viamo	PC-20M
		PC-50M			
	Skin surface transducer	PST-25BT			
		PST-65AT			
		PLT1204ST			
		PST25ST			
		PST50BT			
		PST-30BT			
		PVT-375BT			
		PVT-375SC			
		PVT-382BT			
		PVT712BT			
		PLT-604AT			
PLT-704AT					
PLT-8055AT					
PLT-1202S					
PLT-1204BT					

TRISTEL FUSE FOR STELLA MATERIAL COMPATIBILITY

Manufacturer	Instrument Type	Instrument Series	Model
Toshiba	Skin surface transducer	Aplio & Viamo	PLT-1204BX
			PLT-704SBT
			PLT-705BT
			PLT-1005BT
			PLT704ST
			PLT-1204MV
		Nemio & Famio	PVF-381MT
			PSM-37AT
			PSM-70AT
			PVF-575AT
			PVQ-375A
			PVQ-641V
			PVQ-674A
			PVF-3620ST
			PVM-375AT
			PSM-37CT
			PSM-50AT
			PLG-506M
			PLG-805S
			PLQ-703A
			PLQ-805A
			PLQ-1203A
			PLM-503AT
			PLM-703AT
		PLM-805AT	
		PSM-20CT	
		PSM-30BT	
		Xario	PVU-375BT
			PVU-382BT
			PSU-50BT
PSU-25BT			
PSU-30BT			
PSU-70BT			
PVU712BT			
PLU-704BT			
Unisensor	Manometry catheter	UniTip	Not specified
Xion Medical	Nasopharyngoscopes	Not specified	Not specified

Disclaimer: Material compatibility testing performed via different methods including wiping and soaking. Tested materials subjected to varying number of cycles.

For more information contact Tristel team on cooperation@tristel.com



Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Australian sponsor: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

New Zealand: Tristel New Zealand Limited, P.O Box 15451, Tauranga, New Zealand
T +64 (0)7 5771560 - F +64 (0)7 5771567 - E info@tristel.co.nz

Hong Kong & Taiwan: Tristel Asia Limited, 21st Floor, 168 Electric Road, North Point, Hong Kong
T +852 2895 6968 - F +852 2869 4388

Copyright © Tristel Solutions
Mkt-Val-488-1
May 2018