

All you need is Fuse!



Fuse is all you need: it is the best solution for preventing outbreaks and ward closures

Manufactured by Tristel in the UK Fuse is a high-level disinfectant designed specifically for use on large surface areas such as wall and floors.

Fuse has been extensively tested in accredited laboratories worldwide and is proven effective with a contact time of only five minutes.

No more counting chlorine tablets or calculating dilutions
One sachet of Fuse diluted into five litres of water destroys a wide range of microorganisms including C. diff., MRSA, bloodborne viruses, VRE, CPE, Norovirus, Mycobacterium terrae, MDRAB and E. coli.

Nothing beats Fuse for speed, efficacy, compatibility and ease of use.



40 sachets/box

Visit us at IPS on stand 56 for a free sample!

Already use Fuse?
Enter our on-stand Training Portal cubicle to certify with Tristel or request your unique access code by emailing training@tristel.com

Tristel

T +44 (0)1638 721500

E mail@tristel.com

www.tristel.com

Fuse for Medical Surfaces Fragranced: FAL649
Fuse for Medical Surfaces Unfragranced: FAL650

Use biocides safely. Always read the label, user guide and product information before use.